

Product Datasheet

Mouse Fbg ELISA Kit (orb409472)

Catalog Number	orb409472
Category	Assays and Kits
Description	The Mouse Fbg ELISA Kit detects the target Fbg in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.156 µg/mL-10 µg/mL. The kit offers a sensitivity of 0.039 µg/mL.
Target	Fbg
Reactivity	Mouse
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 µg/mL-10 µg/mL
Sensitivity	0.039 µg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)