

Product Datasheet

Mouse Neutrophil Elastase ELISA Kit (orb409544)

Catalog Number	orb409544
Category	Assays and Kits
Description	The Mouse Neutrophil Elastase ELISA Kit detects the target ELA2 in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 0.312 ng/mL-20 ng/mL. The kit offers a sensitivity of 0.078 ng/mL.
Target	ELA2
Reactivity	Mouse
UniProt ID	Q3UP87
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)