

Product Datasheet

Mouse LTA ELISA Kit (orb409549)

Catalog Number	orb409549
Category	Assays and Kits
Description	The Mouse LTA ELISA Kit detects the target LTA in mouse samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 18.75 pg/mL-1200 pg/mL. The kit offers a sensitivity of 4.7 pg/mL.
Target	LTA
Reactivity	Mouse
UniProt ID	P09225
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	18.75 pg/mL-1200 pg/mL
Sensitivity	4.7 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)