

## Product Datasheet

### Mouse TNFAIP6 ELISA Kit (orb409554)

<b>Catalog Number</b>	orb409554
<b>Category</b>	Assays and Kits
<b>Description</b>	The Mouse TNFAIP6 ELISA Kit detects the target TNFAIP6 in mouse samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 62.5 pg/mL-4000 pg/mL. The kit offers a sensitivity of 15.6 pg/mL.
<b>Target</b>	TNFAIP6
<b>Reactivity</b>	Mouse
<b>UniProt ID</b>	<b>O08859</b>
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Assay Type</b>	Sandwich
<b>Assay Time</b>	1-5h
<b>Range</b>	62.5 pg/mL-4000 pg/mL
<b>Sensitivity</b>	15.6 pg/mL
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)