

Product Datasheet

Pig PYY ELISA Kit (orb409731)

Catalog Number	orb409731
Category	Assays and Kits
Description	The Pig PYY ELISA Kit detects the target PYY in porcine samples. It measures levels in sample types including serum and plasma. The detection range is 62.5 pg/mL-1000 pg/mL. The kit offers a sensitivity of 16.88 pg/mL.
Target	PYY
Reactivity	Porcine
UniProt ID	P68005
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	62.5 pg/mL-1000 pg/mL
Sensitivity	16.88 pg/mL
Sample Types	serum, plasma
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)