

Product Datasheet

Pig Luteinizing hormone ELISA Kit (orb409747)

Catalog Number	orb409747
Category	Assays and Kits
Description	The Pig Luteinizing hormone ELISA Kit detects the target LH in porcine samples. This is a quantitative competitive ELISA kit.
Target	LH
Reactivity	Porcine
UniProt ID	P01232
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039/ml
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)