

Product Datasheet

Pig Immunoglobulin M ELISA Kit (orb409772)

Catalog Number	orb409772
Category	Assays and Kits
Description	The Pig Immunoglobulin M ELISA Kit detects the target IgM in porcine samples. It measures levels in sample types including serum and plasma. The detection range is 0.039 µg/mL-10 µg/mL. The kit offers a sensitivity of 0.039 µg/mL.
Target	IgM
Reactivity	Porcine
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	0.039 µg/mL-10 µg/mL
Sensitivity	0.039 µg/mL
Sample Types	serum, plasma
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)