

Product Datasheet

rabbit FSH ELISA Kit (orb410001)

Catalog Number	orb410001
Category	Assays and Kits
Description	The rabbit FSH ELISA Kit detects the target FSH in rabbit samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 1.56 mIU/mL-100 mIU/mL. The kit offers a sensitivity of 0.39 mIU/mL.
Target	FSH
Reactivity	Rabbit
UniProt ID	Q6IY73
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	1.56 mIU/mL-100 mIU/mL
Sensitivity	0.39 mIU/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)