

Product Datasheet

Rat S100A6 ELISA Kit (orb410467)

Catalog Number	orb410467
Category	Assays and Kits
Description	The Rat S100A6 ELISA Kit detects the target S100A6 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 2.4 ng/mL-150 ng/mL. The kit offers a sensitivity of 0.6 ng/mL.
Target	S100A6
Reactivity	Rat
UniProt ID	P30801
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	2.4 ng/mL-150 ng/mL
Sensitivity	0.6 ng/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)