

Product Datasheet

Rat CD4 ELISA Kit (orb410620)

Catalog Number	orb410620
Category	Assays and Kits
Description	The Rat CD4 ELISA Kit detects the target CD4 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.78 U/mL-50 U/mL. The kit offers a sensitivity of 0.195 U/mL.
Target	CD4
Reactivity	Rat
UniProt ID	P05540
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.78 U/mL-50 U/mL
Sensitivity	0.195 U/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)