

Product Datasheet

Rat S100B ELISA Kit (orb410663)

Catalog Number	orb410663
Category	Assays and Kits
Description	The Rat S100B ELISA Kit detects the target S100B in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 3.12 pg/mL-200 pg/mL. The kit offers a sensitivity of 0.78 pg/mL.
Target	S100B
Reactivity	Rat
UniProt ID	P04631
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)