

Product Datasheet

Rat IL17F ELISA Kit (orb410713)

Catalog Number	orb410713
Category	Assays and Kits
Description	The Rat IL17F ELISA Kit detects the target IL17F in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 15.6 pg/mL-1000 pg/mL. The kit offers a sensitivity of 3.9 pg/mL.
Target	IL17F
Reactivity	Rat
UniProt ID	Q5BJ95
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	15.6 pg/mL-1000 pg/mL
Sensitivity	3.9 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)