

Product Datasheet

Rat MAO ELISA Kit (orb410794)

Catalog Number	orb410794
Category	Assays and Kits
Description	The Rat MAO ELISA Kit detects the target MAO in rat samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 93.75 mU/mL-6000 mU/mL. The kit offers a sensitivity of 23.43 mU/mL.
Target	MAO
Reactivity	Rat
UniProt ID	P21396
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	93.75 mU/mL-6000 mU/mL
Sensitivity	23.43 mU/mL
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)