

## Product Datasheet

### Rat SFTPB ELISA Kit (orb410834)

|                          |                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb410834                                                                                                                                                                                                                                    |
| <b>Category</b>          | Assays and Kits                                                                                                                                                                                                                              |
| <b>Description</b>       | The Rat SFTPB ELISA Kit detects the target SFTPB in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 3.12 pg/mL-200 pg/mL. The kit offers a sensitivity of 0.78 pg/mL. |
| <b>Target</b>            | SFTPB                                                                                                                                                                                                                                        |
| <b>Reactivity</b>        | Rat                                                                                                                                                                                                                                          |
| <b>UniProt ID</b>        | <b>P22355</b>                                                                                                                                                                                                                                |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                                                                                                                                                                                 |
| <b>Assay Type</b>        | Sandwich                                                                                                                                                                                                                                     |
| <b>Assay Time</b>        | 1-5h                                                                                                                                                                                                                                         |
| <b>Range</b>             | 3.12 pg/mL-200 pg/mL                                                                                                                                                                                                                         |
| <b>Sensitivity</b>       | 0.78 pg/mL                                                                                                                                                                                                                                   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                                                                                                                                                                                            |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                        |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)