

Product Datasheet

Rat Galanin ELISA Kit (orb410849)

Catalog Number	orb410849
Category	Assays and Kits
Description	The Rat Galanin ELISA Kit detects the target GAL in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 4.7 pg/mL-300 pg/mL. The kit offers a sensitivity of 1.17 pg/mL.
Target	GAL
Reactivity	Rat
UniProt ID	P10683
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	4.7 pg/mL-300 pg/mL
Sensitivity	1.17 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)