

Product Datasheet

Rat luteinizing hormone ELISA Kit (orb410915)

Catalog Number	orb410915
Description	Rat luteinizing hormone (LH) Elisa kit
Reactivity	Rat
Range	0.3 mIU/ml-60 mIU/ml
Target	LH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Assay Time	1-5h
Uniprot ID	P01230
Sensitivity	0.15 mIU/ml
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)