

Product Datasheet

Rat MMP14 ELISA Kit (orb410933)

Catalog Number	orb410933
Category	Assays and Kits
Description	The Rat MMP14 ELISA Kit detects the target MMP14 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 46.88 pg/mL-3000 pg/mL. The kit offers a sensitivity of 11.75 pg/mL.
Target	MMP14
Reactivity	Rat
UniProt ID	Q10739
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	46.88 pg/mL-3000 pg/mL
Sensitivity	11.75 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)