

Product Datasheet

Rat CSF1 ELISA Kit (orb410973)

Catalog Number	orb410973
Category	Assays and Kits
Description	The Rat CSF1 ELISA Kit detects the target CSF1 in rat samples. It measures levels in sample types including serum, plasma, cell culture supernates and tissue homogenates. The detection range is 0.156 ng/ml - 10 ng/ml. The kit offers a sensitivity of 0.039 ng/ml.
Target	CSF1
Reactivity	Rat
UniProt ID	Q8JZQ0
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.156 ng/ml - 10 ng/ml
Sensitivity	0.039 ng/ml
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)