

Product Datasheet

Rat CC16 ELISA Kit (orb411010)

Catalog Number	orb411010
Category	Assays and Kits
Description	The Rat CC16 ELISA Kit detects the target CC16 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 3.12 pg/mL-200 pg/mL. The kit offers a sensitivity of 0.78 pg/mL.
Target	CC16
Reactivity	Rat
UniProt ID	P17559
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.78 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)