

Product Datasheet

Rat BNP ELISA Kit (orb411047)

Catalog Number	orb411047
Category	Assays and Kits
Description	The Rat BNP ELISA Kit detects the target BNP in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 125 pg/mL-8000 pg/mL. The kit offers a sensitivity of 31.25 pg/mL.
Target	BNP
Reactivity	Rat
UniProt ID	P13205
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	125 pg/mL-8000 pg/mL
Sensitivity	31.25 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)