

Product Datasheet

Rat C10orf59 ELISA Kit (orb411141)

Catalog Number	orb411141
Category	Assays and Kits
Description	The Rat C10orf59 ELISA Kit detects the target C10orf59 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 62.5 ng/mL-4000 ng/mL. The kit offers a sensitivity of 15.6 ng/mL.
Target	C10orf59
Reactivity	Rat
UniProt ID	Q5U2W9
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	62.5 ng/mL-4000 ng/mL
Sensitivity	15.6 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)