

Product Datasheet

Rat RBP4 ELISA Kit (orb411153)

Catalog Number	orb411153
Category	Assays and Kits
Description	The Rat RBP4 ELISA Kit detects the target RBP4 in rat samples. This is a quantitative competitive ELISA kit.
Target	RBP4
Reactivity	Rat
UniProt ID	P04916
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	39 ng/ml-2500 ng/ml
Sensitivity	9.75 ng/ml
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)