

Product Datasheet

Rat Melatonin ELISA Kit (orb411200)

Catalog Number	orb411200
Category	Assays and Kits
Description	The Rat Melatonin ELISA Kit detects the target MT in rat samples. It measures levels in sample types including serum and plasma. The detection range is 12.5 pg/mL-500 pg/mL. The kit offers a sensitivity of 12.5 pg/mL.
Target	MT
Reactivity	Rat
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	12.5 pg/mL-500 pg/mL
Sensitivity	12.5 pg/mL
Sample Types	serum, plasma
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)