

Product Datasheet

Rat ACE ELISA Kit (orb411260)

Catalog Number	orb411260
Category	Assays and Kits
Description	The Rat ACE ELISA Kit detects the target ACE in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 31.25 ng/ml - 2000 ng/ml. The kit offers a sensitivity of 7.81 ng/ml.
Target	ACE
Reactivity	Rat
UniProt ID	P47820
Application notes	Detection Wavelength: 450 nm
Assay Type	Competitive
Assay Time	1-5h
Range	31.25 ng/ml - 2000 ng/ml
Sensitivity	7.81 ng/ml
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)