

Product Datasheet

Rat HTR2A ELISA Kit (orb411282)

Catalog Number	orb411282
Category	Assays and Kits
Description	The Rat HTR2A ELISA Kit detects the target HTR2A in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates and cell lysates. The detection range is 31.25 pg/mL-2000 pg/mL. The kit offers a sensitivity of 7.8 pg/mL.
Target	HTR2A
Reactivity	Rat
UniProt ID	P14842
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	31.25 pg/mL-2000 pg/mL
Sensitivity	7.8 pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)