

Product Datasheet

Rat IGF1R ELISA Kit (orb411303)

Catalog Number	orb411303
Category	Assays and Kits
Description	The Rat IGF1R ELISA Kit detects the target IGF1R in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 3.12 pg/mL-200 pg/mL. The kit offers a sensitivity of 0.8 pg/mL.
Target	IGF1R
Reactivity	Rat
UniProt ID	P24062
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	3.12 pg/mL-200 pg/mL
Sensitivity	0.8 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)