

Product Datasheet

Rat IGFBP2 ELISA Kit (orb411307)

Catalog Number	orb411307
Category	Assays and Kits
Description	The Rat IGFBP2 ELISA Kit detects the target IGFBP2 in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.625 ng/mL-40 ng/mL. The kit offers a sensitivity of 0.156 ng/mL.
Target	IGFBP2
Reactivity	Rat
UniProt ID	P12843
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Sample Types	serum, plasma, tissuehomogenates.
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)