

Product Datasheet

Rat GLP2 ELISA Kit (orb411318)

Catalog Number	orb411318
Category	Assays and Kits
Description	The Rat GLP2 ELISA Kit detects the target GCG in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.78 pg/mL-50 pg/mL. The kit offers a sensitivity of 0.195 pg/mL.
Target	GCG
Reactivity	Rat
UniProt ID	P06883
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.78 pg/mL-50 pg/mL
Sensitivity	0.195 pg/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)