

Product Datasheet

Rat ADH ELISA Kit (orb411323)

Catalog Number	orb411323
Category	Assays and Kits
Description	The Rat ADH ELISA Kit detects the target ADH in rat samples. It measures levels in sample types including serum, plasma and tissue homogenates. The detection range is 0.312 ng/mL-20 ng/mL. The kit offers a sensitivity of 0.078 ng/mL.
Target	ADH
Reactivity	Rat
UniProt ID	P06757
Application notes	Detection Wavelength: 450 nm
Assay Type	Sandwich
Assay Time	1-5h
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Sample Types	serum, plasma, tissue homogenates
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)