

Product Datasheet

DON ELISA Kit (orb411508)

Catalog Number	orb411508
Category	Assays and Kits
Description	The DON ELISA Kit detects the target Deoxynivalenol(DON) in sample types including feed and soybean meal. The detection range is 1 µg/kg~81 µg/kg. The kit offers a sensitivity of 1 µg/kg.
Target	Deoxynivalenol(DON)
Application notes	Detection Wavelength: 450 nm
Assay Time	1-5h
Range	1 µg/kg~81 µg/kg
Sensitivity	1 µg/kg
Sample Types	feed, soybean meal
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)