

Product Datasheet

Mouse CXCL7/PPBP ELISA Kit (orb413112)

Catalog Number	orb413112
Category	Assays and Kits
Description	Mouse CXCL7/PPBP ELISA Kit
Target	Chemokine subfamily B Cys-X-Cys
Reactivity	Mouse
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K40-Y113
UniProt ID	Q9EQI5
Tested applications	ELISA
Specificity	Natural and recombinant mouse CXCL7
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	15.6 pg/ml - 1,000 pg/ml
Sensitivity	2 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Mouse CXCL7 ELISA Kit

Mouse CXCL7 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)