

Product Datasheet

AGXT Antibody (orb414605)

Catalog Number	orb414605
Category	Antibodies
Description	This AGXT Antibody is an unconjugated polyclonal product. It targets AGXT using a Recombinant Human Serine--pyruvate aminotransferase protein (112-294AA) as the immunogen. This antibody is suitable for ELISA, IHC, WB. Purification: >95%, Protein G purified.
Target	AGXT
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Serine--pyruvate aminotransferase protein (112-294AA)
UniProt ID	P21549
Tested applications	ELISA, IHC, WB
Dilution range	WB:1:500-5000, IHC-P:1:20-200

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only



Immunohistochemistry of paraffin-embedded human liver cancer tissue using AGXT antibody



Western blot analysis of A549 whole cell lysate, MCF-7 whole cell lysate using AGXT antibody



Immunohistochemistry of paraffin-embedded human liver tissue using AGXT antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558