

Product Datasheet

HIST1H3A (Ab-27) Antibody (orb416465)

Catalog Number	orb416465
Category	Antibodies
Description	This HIST1H3A (Ab-27) Antibody is an unconjugated polyclonal product. It targets HIST1H3A using a Peptide sequence around site of Lys (27) derived from Human Histone H3.1 as the immunogen. This antibody is suitable for ELISA, IF, IHC. Purification: Antigen Affinity Purified.
Target	HIST1H3A
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	Antigen Affinity Purified
Immunogen	Peptide sequence around site of Lys (27) derived from Human Histone H3.1
UniProt ID	P68431
Tested applications	ELISA, IF, IHC
Dilution range	IHC-P:1:20-200, IF:1:50-200

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only



Immunohistochemistry of paraffin-embedded human lung cancer tissue using HIST1H3A (Ab-27) antibody



Immunofluorescent analysis of HeLa cells using HIST1H3A (Ab-27) antibody



Immunohistochemistry of paraffin-embedded human colon cancer tissue using HIST1H3A (Ab-27) antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558