

Product Datasheet

Mono/Di/Tri-methyl-Histone H3.1 K79 Recombinant Monoclonal Antibody (orb418656)

Description	Rabbit monoclonal antibody to Mono/Di/Tri-methyl-Histone H3.1 (K79)
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A synthesized peptide
Target	HIST1H3A
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Clone Number	2G10
Antibody Type	Recombinant Antibody
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)