

Product Datasheet

Human TICAM1 protein (orb419267)

Description	Recombinant Human TIR domain-containing adapter molecule 1
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	Tag-Free
Note	For research use only
Application notes	Tag Info: NO-taggedExpression Region: 1-712aaSequence Info: Full Length
Protein Sequence	MACTGPSLPSAFDILGAAGQDKLLYLKHKLKTPRPGCQGQDLLHAMVLLKLGQETEARISLE ALKADAVARLVARQWAGVDSTEDPEEPPDVSWAVARLYHLLAEKLCPASLRDVAYQEAV RTLSSRDDHRLGELQDEARNRCGWDIAGDPGSIRTLQSNLGCLPPSSALPSGTRSLPRPID GVSDWSQGCSLRSTGSPASLASNLEISQSPTMPFLSLHRSPHGPSKLCDDPQASLVPEPVP GGCQEPEEMSWPPSGEIASPPELPSSPPPGLEVPAPDATSTGLPDTAAPETSTNYPVECTE GSAGPQSLPLPILEPVKNPCSVKDQTPQLQSVEDTTSPNTKPCPTPTTPETSPPPPPPPPSST PCSAHLTPSSLPSSLESSSEQKFYNFVILHARADEHIALRVREKLEALGVPDGFATCEDFQV PGRGELSCLQDAIDHSAFIILLTSNFDCLSLHQVNQAMMSNLTRQGSPDCVIPFLPLESSP AQLSSDTASLLSGLVRLDEHSQIFARKVANTFKPHRLQARKAMWRKEQDTRALREQSQHL DGERMQAAALNAAYSAYLQSYLSYQAQMEQLQVAFGSHMSFGTGAPYGARMPFGGQVPL GAPPPFPTWPGCPQPPPLHAWQAGTPPPSPQPAAFPQSLPFPQSPAFPTASPAPPQSPGL QPLIIHHAQMVLQGLNNHMMWNQRGSQAPEDKTQEAE
Purity	Greater than 90% as determined by SDS-PAGE.


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Source	in vitro E.coli expression system
MW	76.4 kDa
Uniprot ID	Q8IUC6
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)