

Product Datasheet

CD11b Antibody, FITC (orb423282)

Description	Mouse monoclonal antibody against CD11b conjugated to FITC
Reactivity	Human
Conjugation	FITC
Immunogen	Purified human PBL Monocytes
Preservatives	1mg/ml in PBS (after reconstitution).
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: Blocking and staining antibody. For staining, use 10µl/1,000,000 cells. Titer for blocking LPS binding should be determined by the investigator
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	hCD11b
Solubility (25°C)	Reconstitute with of H2O. Mix gently, wash the sides of the vial and wait 30-60 seconds before use.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)