

Product Datasheet

Human PMP2 (His) Protein (orb427753)

Catalog Number	orb427753
Category	Proteins
Description	Recombinant of human PMP2 (His) protein
Form/Appearance	Sterile Filtered colorless liquid formulation.
Buffer/Preservatives	PMP2 His-Tag is supplied in 20mM Tris HCl pH-8 and 50% glycerol.
Purity	Greater than 95% as determined by SDS-PAGE.
Application notes	Recombinant & Natural Proteins
Source	Escherichia Coli
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)