

Product Datasheet

DYKDDDDK TAG Antibody (orb430674)

Catalog Number	orb430674
Category	Antibodies
Description	Rabbit polyclonal antibody to DYKDDDDK TAG
Target	DYKDDDDK TAG
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	Preservative: 0.1% Sodium Azide (NaN ₃); Buffer: Phosphate buffered saline
Purity	Purified
Immunogen	Keyhole limpet haemocyanin conjugated Enterokinase cleavage site peptide DYKDDDDK. Residues of cysteine and glycine were used to facilitate coupling at the C-terminal end.
Tested applications	ELISA, WB
Specificity	DYKDDDDK TAG
Antibody Type	Secondary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)