

Product Datasheet

SIGIRR Antibody (orb431038)

| | |
|-----------------------------|--|
| Catalog Number | orb431038 |
| Category | Antibodies |
| Description | Rabbit polyclonal antibody to SIGIRR |
| Target | SIGIRR |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | Polyclonal IgG |
| Conjugation | Unconjugated |
| Reactivity | Mouse |
| Form/Appearance | Purified IgG - liquid |
| Concentration | 1.0 mg/ml |
| Buffer/Preservatives | Preservative: 0.02% Sodium Azide (NaN ₃); Buffer: Phosphate buffered saline |
| Purity | Purified |
| Immunogen | A synthetic peptide corresponding to 15 amino acids near the carboxy terminus of mouse SIGIRR. |
| UniProt ID | Q6IA17, Q9JLZ8 |
| Tested applications | IH-Fr, WB |
| Application notes | (C-TERMINAL) |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|-------------------------|---|
| Specificity | SIGIRR |
| Cross Reactivity | Human,Rat |
| Antibody Type | Primary Antibody |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)