

Product Datasheet

GLUTAMATE DECARBOXYLASE 1/2 Antibody (orb431732)

Catalog Number	orb431732
Category	Antibodies
Description	Rabbit polyclonal antibody to GLUTAMATE DECARBOXYLASE 1/2
Target	GLUTAMATE DECARBOXYLASE 1/2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Reactivity	Rat
Form/Appearance	Purified IgG - liquid
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃); Buffer: Phosphate buffered saline
Purity	Purified
Immunogen	Synthetic peptide sequence D-F-L-I-E-E-I-E-R-L-G-Q-D-L from the C terminus of rat GAD2
UniProt ID	Q05683
Tested applications	ELISA, IF, IH-Fr, IHC-P, WB
Application notes	65/67kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Specificity	GLUTAMATE DECARBOXYLASE 1/2
Cross Reactivity	Mouse,Human,Cat
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)