

Product Datasheet

RABBIT ANTI SHEEP INTERLEUKIN-1 BETA (orb431755)

Catalog Number	orb431755
Description	Rabbit polyclonal antibody to IL-1 BETA
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant IL-1 beta.
Target	IL-1 BETA
Preservatives	Serum - liquid
Form/Appearance	Serum - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Serum
Uniprot ID	P21621
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)