

Product Datasheet

CD184/CXCR4 Antibody (orb431766)

Catalog Number	orb431766
Category	Antibodies
Description	Rabbit polyclonal antibody to CD184 / CXCR4
Target	CD184 / CXCR4
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Buffer/Preservatives	Preservative: 0.02% Sodium Azide (NaN ₃); Buffer: Phosphate buffered saline
Purity	Purified
Immunogen	Peptide sequence corresponding to amino acids 1-14 of human CD184 / CXCR4.
UniProt ID	P61073
Tested applications	FC, IH-Fr, IP, WB
Application notes	(N-TERMINAL)
Specificity	CD184 / CXCR4

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)