

Product Datasheet

RABBIT ANTI BOVINE INTERLEUKIN-1 BETA:Biotin (orb431950)

Description	Rabbit polyclonal antibody against IL-1 BETA conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant bovine IL-1 beta.
Target	IL-1 BETA
Preservatives	Purified IgG conjugated to Biotin - lyophilized
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.25mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P09428
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)