

Product Datasheet

V5-TAG Antibody (Biotin) (orb432441)

Catalog Number	orb432441
Category	Antibodies
Description	Mouse monoclonal antibody against V5-TAG conjugated to Biotin
Target	V5-TAG
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a
Conjugation	Biotin
Reactivity	Virus
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	1.0 mg/ml
Buffer/Preservatives	Preservative: 0.09% sodium azide (NaN ₃); Buffer: Phosphate buffered saline
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
UniProt ID	P11207
Tested applications	ELISA, IH-Fr, WB
Specificity	V5-TAG
Antibody Type	Secondary Antibody
Clone Number	SV5-Pk1

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)