

## Product Datasheet

### MOUSE ANTI HUMAN HLA G:FITC (orb433118)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against HLA G conjugated to FITC                     |
| <b>Species/Host</b>        | Mouse                                                                          |
| <b>Reactivity</b>          | Human                                                                          |
| <b>Conjugation</b>         | FITC                                                                           |
| <b>Tested Applications</b> | FC                                                                             |
| <b>Immunogen</b>           | Recombinant human HLA-G refolded with beta 2 microglobulin.                    |
| <b>Target</b>              | HLA G                                                                          |
| <b>Preservatives</b>       | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| <b>Concentration</b>       | 0.1 mg/ml                                                                      |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                       |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | IgG1                                                                           |
| <b>Clonality</b>           | Monoclonal                                                                     |
| <b>Clone Number</b>        | MEM-G/9                                                                        |
| <b>Uniprot ID</b>          | <b>P17693</b>                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)