

Product Datasheet

CD284 Antibody (FITC) (orb433149)

| | |
|-----------------------------|----------------------------------------------------------------------------------------------------------------|
| Catalog Number | orb433149 |
| Category | Antibodies |
| Description | Mouse monoclonal antibody against CD284 conjugated to FITC |
| Target | CD284 |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Isotype | IgG2a |
| Conjugation | FITC |
| Reactivity | Human |
| Form/Appearance | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| Concentration | 0.1 mg/ml |
| Buffer/Preservatives | Preservative: 0.09% sodium azide (NaN ₃)1% bovine serum albumin; Buffer: Phosphate buffered saline |
| Immunogen | Ba/F3 cell line expressing TLR4 (CD284). |
| UniProt ID | O00206 |
| Tested applications | FC |
| Specificity | CD284 |
| Cross Reactivity | Bovine,Dog,Guinea Pig,Pig,Rhesus Monkey |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antibody Type | Primary Antibody |
| Clone Number | HTA125 |
| Storage | This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended. This product is photosensitive and should be protected from light. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)