

Product Datasheet

MOUSE ANTI SHEEP CD4:RPE (orb433338)

Catalog Number	orb433338
Category	Antibodies
Description	Mouse monoclonal antibody against CD4 conjugated to RPE
Target	CD4
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a
Conjugation	RPE
Reactivity	Sheep
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Immunogen	Fetal thymocytes.
UniProt ID	P05542
Tested applications	FC
Cross Reactivity	Goat
Clone Number	44.38
Storage	+4°C DO NOT FREEZE
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)