

Product Datasheet

CD45RO Antibody (RPE) (orb433697)

Catalog Number	orb433697
Category	Antibodies
Description	Mouse monoclonal antibody against CD45RO conjugated to RPE
Target	CD45RO
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG3
Conjugation	RPE
Reactivity	Bovine
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin5% Sucrose; Buffer: Phosphate buffered saline
Immunogen	Bovine peripheral blood mononuclear cells.
Tested applications	FC
Specificity	CD45RO
Cross Reactivity	Goat
Antibody Type	Primary Antibody
Clone Number	IL-A116

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Note

For research use only

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)