

Product Datasheet

MOUSE ANTI RABBIT LIGHT CHAIN:HRP (orb434963)

Catalog Number	orb434963
Category	Antibodies
Description	Mouse monoclonal antibody against LIGHT CHAIN conjugated to HRP
Target	LIGHT CHAIN
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2b
Conjugation	HRP
Reactivity	Rabbit
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Buffer/Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Immunogen	Rabbit light chain
Tested applications	ELISA, WB
Clone Number	SB62a
Storage	-20°C only
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)