

Product Datasheet

MOUSE ANTI SHEEP CD230 (aa143-154) (orb434964)

Catalog Number	orb434964
Description	Mouse monoclonal antibody to CD230
Species/Host	Mouse
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	A truncated form of recombinant Ovine PRP spanning residues aa94 - aa233.
Target	CD230
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	(aa143 - aa154)
Isotype	IgG1
Clonality	Monoclonal
Clone Number	ROS-BH1

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Purity	Purified
Uniprot ID	P10279, P23907
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)